




SPECIFICATION SHEET

3277DK

SUPER GRIP DISPOSABLE GLOVES

FIT, DEXTERITY, TRACTION GRIP, AMBIDEXTROUS

- Embossed diamond pattern on the outside of the glove for ergonomic traction and exceptional wet grip.
- Resistant to Fentanyl
- Tested for use with Chemotherapy drugs
- Medical Grade
- Chemical resistant 100% 8 mil nitrile
- Beaded cuff for secure fit around the arm
- CE level 5 dexterity for superior fit, feel, and comfort
- Powder Free
- 9.5" Length
- Ambidextrous fit for less waste
- Industries: Industrial, Construction, Automotive, Agriculture



Technical Data

Material	100 % Nitrile
Standards	<p>ASTM D6978-05 - Resistance to Permeation of Permeation by Chemotherapy Drugs and Fentanyl</p> <p>ISO 374-1:2016 Type B - Chemical Resistance</p> <p>EN 374-2:2014 - Penetration Resistance</p> <p>EN 16523-1:2015 - Permeation Resistance</p> <p>EN 374-4:2013 - Resistance to Degradation by Chemicals</p> <p>EN ISO 374-5:2016 - Micro-organism Protection</p> <p>EN 455-1:2000 - Medical Glove for Single Use</p> <p>FDA CFR 21 177.2600 - Suitable for Food Contact</p>
Cuff Closure	Beaded Cuff
Mil Thickness	8 Mil
Colors	Black
Sizes	S-XXL
Packaging	100 gloves/box 10 boxes/carton

Embossed diamond pattern provides superior grip



Ordering Information

Item	Siz	UPC
3277DK/ 8	e S	73179117390
3277DK/ 9	M	0
3277DK/10	L	73179117391
3277DK/11	X1	7
3277DK/12	X2	73179117392
		4
		73179117393
		1
		73179117394
		8



All safety equipment should be selected in conjunction with an approved safety program that assesses your specific hazards. The selection and use of specific products sold by Majestic is exclusively within the control of, and remains the sole responsibility of the customer. Many variables beyond the control of Majestic affect the results obtained in applying these products. Inappropriate selection and/or use of safety products may cause serious injury.



P: 724-779-2278

www.sq1med.com